




Stainless Steel Impression Trays Made of 18/8 Stainless Steel

**Solid Partial Set of 10
225-001S**



786-1084
U1

786-1085
U2

786-1086
U3

786-1087
U4

786-1088
U5

786-1084
L1


786-1085
L2

786-1086
L3

786-1087
L4

786-1088
L5

All the above impression trays are available in these 5 sizes





UPPER	UPPER				
	U1	U2	U3	U4	U5
A	62	60	53	54	52
B	83	75	73	69	68

All Dimension are in mm

LOWER	LOWER				
	L1	L2	L3	L4	L5
A	63	61	57	60	53
B	82	81	78	74	71

All Dimension are in mm



786-1089
L1

786-1090
L2

786-1091
L3

786-1092
L4

786-1093
L5

786-1089
U1

786-1090
U2

786-1091
U3

22786-1092
5-009
U4


786-1093
U5

**Perforated Partial Set of 10
225-010S**




Stainless Steel Impression Trays Made of 18/8 Stainless Steel

**Solid Full Denture Set of 8
225-011S**



All the above impression trays are available in these 4 sizes




UPPER		XL	L	M	S
	A	70	67	64	60
B	79	76	70	73	

All Dimension are in mm

LOWER		XL	L	M	S
	A	64	64	60	57
B	86	70	73	70	

All Dimension are in mm



**Perforated Full Denture Set of 8
225-018S**




Stainless Steel Impression Trays Made of 18/8 Stainless Steel


**Solid Set of 6
225-019S**




786-1102
U1




786-1103
U2




786-1104
U3



786-1102
L1




786-1103
L2



786-1104
L3

All the above impression trays are available in these 3 sizes



UPPER


	U1	U2	U3
A	67	64	57
B	81	71	70

All Dimension are in mm


LOWER

	L1	L2	L3
A	60	57	54
B	81	79	76

All Dimension are in mm



786-1105
U1



786-1106
U2



786-1107
U3



786-1105
L1



786-1106
L2




786-1107
L3


**Perforated Set of 6 L2
225-024S**



Stainless Steel Impression Trays Made of 18/8 Stainless Steel



All the above impression trays are available in these 4 sizes




UPPER	Size	A	B	C
	0	67	50	12
	1	72	51	13
	2	78	55	15
	3	80	58	15

All Dimension are in mm

LOWER	Size	A	B	C
	0	70	50	14
	1	77	51	15
	2	79	55	15
	3	80	60	17

All Dimension are in mm





Stainless Steel Impression Trays Made of 18/8 Stainless Steel

**Solid Full Denture Set of 8
225-028S**



786-1111
XL



786-1112
L



786-1113
M



786-1114
S



786-1111
XL



786-1112
L




786-1113
M



786-1114
S

All the above impression trays are available in these 4 sizes




		XL	L	M	S
UPPER	A	70	67	64	60
	B	79	76	70	73

All Dimension are in mm

		XL	L	M	S
LOWER	A	64	64	60	57
	B	86	70	73	70

All Dimension are in mm



786-1115
XL




786-1116
L



786-1117
M



786-1118
S




786-1115
XL



786-1116
L



786-1117
M



786-1118
S


**Perforated Full Denture Set of 8
225-035S**




Stainless Steel Impression Trays Made of 18/8 Stainless Steel




786-1119
U0



786-1120
U1




786-1121
U2



786-1122
U3



786-1119
L0



786-1120
L1



786-1121
L2




786-1122
L3

**Upper Sizes 0,1,2,3,
Depth of the tray 15,9 mm
Height of palate 12,5 mm**


**Lower Sizes 0,1,2,3,
Depth of the tray 12,9 mm**



786-1123
U0



786-1124
U1



786-1125
U2



786-1126
U3



786-1123
L0



786-1124
L1



786-1125
L2




786-1126
L3

**Perforated Full Denture Set of 8
225-043S**



Stainless Steel Impression Trays Made of 18/8 Stainless Steel

Impression trays for crown & bridge work



786-1127
U1P



Solid

786-1127
L1P




786-1128
U1P




Perforated

786-1128
L1P




786-1129
For upper right jaw
For lower left jaw

Perforated



786-1130
For upper left jaw
For lower right jaw



786-1131
Adjustable
Perforated



786-1131A
Adjustable
Solid

Available in 4 Sizes

- 1..... 44 mm
- 2.....55 mm
- 3.....60 mm
- 4.....65 mm

Also available in Solid pattern



Stainless Steel Instrument Trays

Made of 18/8 Stainless Steel


786-1132S

Complete instrument tray with perforated base, lid & insert frame for 16 instrument



786-1132

Lid; dimensions mm, 288x187x29



786-1133

Perforated base; dimensions 284x183x17 mm


786-1125S

Complete instrument tray with solid base, lid & insert frame for 16 instrument



786-1134

Lid; dimensions 288x187x29 mm



786-1135

Solid base; dimensions 284x183x17 mm


786-1137S

Complete instrument tray with perforated base, lid & insert frame for 23 instrument



786-1136

Lid; dimensions mm, 288x187x29



786-1137

Perforated base; dimensions 284x183x17 mm


786-1139S

Complete instrument tray with solid base, lid & insert frame for 23 instrument



786-1138

Lid; dimensions 288x187x29 mm




786-1139

Solid base; dimensions 284x183x17 mm




Stainless Steel Impression Trays Made of 18/8 Stainless Steel

Scalers Tray



786-1140




786-1141
6" x 10"
without lid




786-1142
6" x 8"
with lid




786-1143
6" x 8"
8" x 10"
10" x 12"
with lid



786-1144




786-1145
Endodontic box
Also available in Solid pattern



786-1146
8" x 12"
with lid




Stainless Steel Instrument Trays Made of 18/8 Stainless Steel



786-1147

Deep

2.25" x 2"
4.75" x 2"
5.5" x 2.5"
6.25" x 2.75"
7" x 3.5"



Kidney trays

786-1148
12"-10"- 8"- 6"




786-1149
6" x 9" - 8" x 12"




Shallow

786-1150
10" x 10"



786-1151
Cotton holder



Round bowls

80x35 mm
125x60 mm

786-1152



786-1153
Spirit lamp



786-1154
Spirit lamp



786-1155
Forcep Jar

25 x 80 mm
30 x 100 mm
100 x 200 mm
75 x 173 mm