








Tooth Extracting Forceps American Pattern




786-192
Fig. 1
Upper centrals, laterals & cuspids



786-193
Fig. 5
Lower third molars



786-194
Fig. 6
Splitting forceps



786-195
Fig. 10 H
upper molars
hook - handle




786-196
Fig. 10 S
upper molars
straight - handle

786-197
Fig. 15
Lower molars



Tooth Extracting Forceps American Pattern



786-199
Fig. 16
Lower molars
(cow-horn beak)

786-201
Fig. 18L
Upper molars
Left side

786-200
Fig. 17
Lower molars


786-198
Fig. 16 S
Lower molars
Children

786-202
Fig. 18R
Upper molars
Right side


786-203
Fig. 23
Lower molars
(Cow-horn beak)




Tooth Extracting Forceps American Pattern




786-204
Fig. 32
Upper molars
& pre molars
strong beak




786-205
Fig. 53 L
Upper molars
Left side



786-206
Fig. 53 R
Upper molars
Right side



786-207
Fig. 62
Universal forceps
for upper teeth & roots




786-208
Fig. 65
upper incisors & root




786-209
Fig. 69
Fragments & small roots
Upper & Lower




Tooth Extracting Forceps American Pattern




786-210
Fig. 88R
Right upper molars



786-211
Fig. 88L
Left upper molars




786-212
Fig. 99C
upper premolars,
molar & incisors




786-213
Fig. 101
Universal forceps
for premolars.



786-214
Fig. 150
Upper premolars,
incisors & roots




786-215
Fig. 150A
Upper premolars,
incisors & roots




786-216
Fig. 150S
Upper premolars,
incisors & roots
for children




Tooth Extracting Forceps American Pattern




786-217
Fig. 151
Lower premolars, incisors & roots




786-218
Fig. 151A
Lower premolars, incisors & roots




786-219
Fig. 151S
Lower premolars, incisors & roots
for children




786-220
Fig. 210H
Upper wisdoms
with hook handle



786-221
Fig. 210S
Upper wisdoms
with straight handle



786-222
Fig. 222
Lower wisdoms



786-223
Fig. 286
Upper roots, premolars & incisors